

 CORIELL INSTITUTE <small>FOR MEDICAL RESEARCH</small>	DNA QUANTIFICATION INFOR FOR SHIPPED DNA	Implementation Date: 5/22/24
	WI 0700-05	Revision Number: I-052224

Please see our website (catalog.coriell.org) for further information.

1. How is DNA extracted?

DNA is extracted either by a Miller’s Salting out procedure or by a magnetic bead purification method. If you require the specific extraction method, please inquire through Coriell Customer Service (customerservice@coriell.org).

2. How is DNA quality control performed?

- DNA concentration is determined using either spectrophotometric absorbance at A260 using the convention that 1 O.D. is equivalent to 50 µg/ml of double-stranded DNA or the Qubit dsDNA BR assay, a dye-based fluorescence method, depending on the repository. Please see #9 for specific information. Samples are read at least two times to verify reading. Instruments are tested quarterly and calibrated as necessary.
- To assess DNA integrity, DNA is evaluated by the Agilent TapeStation. The TapeStation software determines a DNA Integrity Number (DIN) as a measure of gDNA integrity. (Note: This is not performed for all repositories.)
- DNA sample identity is confirmed using a multiplex PCR assay for 6 autosomal microsatellite markers. Gender is determined in the same reaction using an additional primer pair designed to amplify a region of allelic difference between the X and Y–chromosome amelogenin genes. Details of this assay and its importance in the quality control process are discussed below. This assay provides several quality control parameters as follows:
 - It qualifies the DNA as a substrate for PCR.
 - It provides a DNA fingerprint that is specific for the sample and the subject from which it is derived.
 - It determines the gender of the subject from which the DNA is derived.
- The criteria for acceptance of the DNA sample based upon all quality control procedures are as follows:
 - Samples have a 260/280 ratio greater than or equal to 1.65 and less than 2.1.
 - Confirmed identity of sample donor.
 - DNA must be intact, be of high molecular weight, and have a DIN score ≥ 7.0 as determined by the Agilent TapeStation (if specified by repository).

3. What is DNA dissolved in?

DNA is dissolved in TE buffer (10 mM TRIS pH 8.0, 1 mM EDTA, pH 8.0).

4. How much DNA is shipped in the vial?

The amount of DNA shipped in the vial is determined by the repository, and the method of quantitation also differs by repository. Please see #9 for specific information. The quantity of DNA is included on the catalog page for each item as well as the pricing page for each collection.

 CORIELL INSTITUTE <small>FOR MEDICAL RESEARCH</small>	DNA QUANTIFICATION INFO FOR SHIPPED DNA	Implementation Date: 5/22/24
	WI 0700-05	Revision Number: I-052224

5. What is the concentration and volume of the DNA sample?

Most samples shipped by Coriell are between 300 – 375 ng/μl, but this differs by repository. The concentration for each DNA sample is available on the Concentration Sheet that accompanies all shipments. Please note that concentration is determined by different methods for different repositories. Please see #9 for specific information.

The volume in each DNA sample is determined by the following formula: Volume= Quantity of DNA/concentration

7. How is DNA shipped?

DNA in tubes is shipped at ambient temperature and is stable. Please freeze at -20°C upon arrival if the sample will not be used shortly. DNA plates are shipped on dry ice.

8. How should DNA be stored?

DNA samples can be stored indefinitely at -80°C. Repeated freeze/thaw cycles should be avoided to reduce sample degradation. We recommend that you aliquot into sample sizes that will be adequate for future use in the lab.

9. How should I retest the DNA concentration?

We recommend using a Nanodrop or conventional cuvette-based spectrophotometer with TE as the appropriate sample blank. DNA concentration can be calculated from the A260 value using the convention that 1 O.D. is equivalent to 50 μg/ml of double-stranded DNA.

Qubit® dsDNA BR Assay (Q32850; ThermoFisher Scientific) can also be used to determine DNA concentration. This assay uses an ultra-sensitive fluorescent nucleic acid stain for quantifying double-stranded DNA (dsDNA) in solution with a standard fluorometer and fluorescein excitation and emission wavelengths.

Coriell DNA samples have been quantitated either by A260 or, after August 1, 2011, by the Qubit assay (refer to the table below). Please see the Concentration sheet included with your DNA to find out the date of extraction.

If you require that concentration be determined by an additional method by Coriell (or require a Certificate of Analysis), please inquire through Coriell Customer Service (customerservice@coriell.org).

Repository	Method of DNA Quantitation
NINDS	A260
NIGMS	Qubit as of August 1, 2011
NHGRI	Qubit as of August 1, 2011
NIA	Qubit as of August 1, 2011
NEI-AREDS	A260
CDC, HuRef	Qubit as of August 1, 2011

10. How are the DNA vials labeled?

Most vials will have a barcode label with lot number and the date extracted.