 CORIELL INSTITUTE <small>FOR MEDICAL RESEARCH</small>	Human Induced Pluripotent Stem Cell (iPSC) Handling Protocols: Matrigel and mTeSR or other medium	Implementation Date: 9/11/24
	WI 1300-05	Revision Number: J-091124

General Guidelines for Handling Human iPSC cells

- Please refer to the cell line specific certificate of analysis for specific culture instructions including growth media, dissociation reagent, and passaging recommendations. This document is a general guideline for reference only.
- iPSC are cryopreserved in plastic cryovials and shipped on dry ice. If storing the iPSC before thawing, store in liquid nitrogen vapor. Storage directly in liquid nitrogen may result in cracking of o-rings.
- It is highly recommended that a small number of vials are cryopreserved as a master stock before beginning any experimentation.

Media and Reagents

Matrigel® (hESC-qualified, Corning 354277)

Rho Kinase (ROCK) Inhibitor (Y-27632)

Table 1. Media

human iPSC growth medium OPTIONS	human iPSC cryopreservation medium
mTeSR™1	90% Knock-Out Serum Replacement
mTeSR™ Plus	10% DMSO
StemFlex	


	Human Induced Pluripotent Stem Cell (iPSC) Handling Protocols: Matrigel and mTeSR or other medium	Implementation Date: 9/11/24
	WI 1300-05	Revision Number: J-091124

Table 2. Dissociation Reagent – Reagents for dissociation should be chosen based on downstream assay or **information supplied on the Certificate of Analysis** supplied with each iPSC line.

Dissociation Reagent	Concentration	Incubation Period	Stop Medium
Versene (1X)	N/A	At room temperature for 2 minutes or until edges of colonies begin to roll up.	Growth medium
ReLeSR™	N/A	See manufacturer instructions	Growth medium


Coating Plates with Matrigel

1. Reconstitute aliquot of Matrigel with cold medium and mix well. **(This should be performed according to manufacturer’s instructions; there will be lot to lot variation)**
2. Immediately use Matrigel solution to coat 6-well plates.
3. Gently rock plate to spread the Matrigel solution evenly across the plate.
4. Plates should remain properly stored at room temperature for 1 hour to use immediately or seal with parafilm and store at 2-8°C for up to 7 days. To use plates after storage at 2-8°C, remove parafilm and incubate for 1 hour at room temperature.
5. Aspirate Matrigel and add warm mTeSR1 or media specified on the CofA (growth medium).

Thawing Human iPS cells

Perform the following steps in a time-efficient manner to obtain optimal cell viability.

1. Remove iPSC from liquid nitrogen vapor or dry ice and immerse the cryovial in a 37°C water bath. Thaw quickly by gently swirling until only a small piece of frozen material remains. Spray the vial with 70% ethanol before transferring to a biological safety cabinet.
2. Gently add the thawed cell suspension dropwise to 10 ml warm growth medium in a sterile 15 ml conical tube and gently mix cells by swirling.
3. Centrifuge conical tube containing cells at 1100 rpm for 3 min at room temperature.
4. Aspirate supernatant and gently re-suspend cells into 2 ml of warm growth medium

	Human Induced Pluripotent Stem Cell (iPSC) Handling Protocols: Matrigel and mTeSR or other medium	Implementation Date: 9/11/24
	WI 1300-05	Revision Number: J-091124

supplemented with 10 μ M Y-27632 Rho kinase (ROCK) inhibitor (Y-compound).


5. Plate cells in 1 well of 6-well plate (unless stated otherwise on the cell line certificate of analysis) coated with Matrigel and gently rock plate to evenly distribute cells.
6. Replace with fresh medium without Y-compound the day after thaw. It is important to remove the Y-compound from the medium within 20 hours.
7. Maintain cells by daily medium exchange. **WARNING:** Failure to replace medium daily can result in spontaneous differentiation.
8. Colonies should be observed within 1-3 days.

Passaging of Human iPSC

Passage iPSCs at ~75-85% confluency or when colonies approach borders of an adjacent colony. Ideally, iPSCs should be passaged before individual colonies begin differentiating in the center of colony. To avoid spontaneous differentiation, do not allow colonies to overgrow.

Split ratios range from 1:3 to 1:6 or higher and **are indicated on the Certificate of Analysis** for each iPSC line. We recommend that you record the passage information for each cell line to determine the exact growth kinetics of a particular line in your laboratory.

1. Remove spent medium from culture and rinse with PBS.
2. Add 1 ml/well of dissociation reagent for specified amount of time or until edges of colonies begin to roll up (see Table 2 and refer to the CofA for each cell line)
3. Aspirate dissociation reagent.
4. Wash each well with 1ml of growth medium (see Table 2).
5. Add 1 ml of warm growth medium to cells and, using a cell scraper, gently dislodge cells from plate. We recommend working with a 5ml or a 2ml pipette but NOT micropipettes since the small bore of pipette tips can break colony clumps to smaller size.
6. Transfer the detached cell aggregates to a 15 ml conical tube. Rinse each well with 1 ml


	Human Induced Pluripotent Stem Cell (iPSC) Handling Protocols: Matrigel and mTeSR or other medium	Implementation Date: 9/11/24
	WI 1300-05	Revision Number: J-091124

of growth medium to collect any remaining aggregates and add to conical tube containing cells.

7. Re-suspend cells in appropriate volume of growth medium for the appropriate cell density for cell culture vessel. We recommend passaging the cells without centrifugation ('no spin' method).
8. Seed cells onto prepared Matrigel coated plates. Rock plate gently to achieve a uniform cell distribution.

Cryopreservation of Human iPSC


1. Remove spent medium from culture and rinse with PBS.
2. Add 1 ml/well of dissociation reagent for specified amount of time or until edges of colonies begin to roll up (see Table 2).
3. Aspirate dissociation reagent.
4. Wash each well with 1ml of specified stop medium (see Table 2).
5. Add 1 ml of warm growth medium to cells and, using cell scraper, gently dislodge cells from plate.
6. Transfer the detached cell aggregates to a conical tube. Rinse each well with 1 ml of growth medium to collect any remaining aggregates and add to conical tube containing cells.
7. Centrifuge conical tube containing cells at 1100 rpm for 3 min at room temperature.
8. Remove supernatant and re-suspend cells in pre-cooled cryopreservation medium from Table 1 (1 ml for each well of a 6 well plate to be harvested).
9. Transfer 1 ml of cell suspension to each cryovial on ice.
10. Place cryovials in isopropanol freezing container and store at -80°C overnight.

 CORIELL INSTITUTE <small>FOR MEDICAL RESEARCH</small>	Human Induced Pluripotent Stem Cell (iPSC) Handling Protocols: Matrigel and mTeSR or other medium	Implementation Date: 9/11/24
	WI 1300-05	Revision Number: J-091124

11. Transfer vials to liquid nitrogen vapor.

Troubleshooting Tips

Problem	Observation of problem	Possible causes
Spontaneous differentiation	Morphology of differentiated cells can vary but is commonly characterized by <ol style="list-style-type: none"> 1. hypertrophic colonies 2. colonies without distinct borders 3. flattened cells 	<ol style="list-style-type: none"> 1. Low confluency or suboptimal passaging of cells (see "low viability after passage") 2. Poor or inappropriate Matrigel concentration/quality- this can be prevented by testing Matrigel prior to using in an
Non-uniform distribution of colonies within culture vessel	Areas within culture vessel with highly confluent iPSC colonies AND areas with few or no iPSC colonies	Failure to properly rock plates after plating cells. Rock plates back and forth gently immediately following plating - usually rocking plates back and forth and then side to side produces a fairly uniform distribution
Low Viability after Passage	Little to no cell colonies are visible 24 hours after passage	During passaging, clusters may have been disrupted into a single cell suspension-decrease pipetting of cells during passaging and decrease incubation time with Versene.

	Human Induced Pluripotent Stem Cell (iPSC) Handling Protocols: Matrigel and mTeSR or other medium	Implementation Date: 9/11/24
	WI 1300-05	Revision Number: J-091124

Low Viability after Recovery from Cryopreservation	Little to no cell colonies are visible within 4 days after recovery	Lack of or insufficient Rho kinase (ROCK) inhibitor in thawing medium.
--	---	--




Image 1: Typical image of a human induced pluripotent stem cell (iPSC) colony. Note the distinct borders of the colony and morphology of the iPSCs.




Image 2: Image of a differentiated colony. Note the flattened morphology of the cells and lack of distinct borders.